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```
```sql
CREATE DATABASE lims_db CHARACTER SET utf8mb4 COLLATE utf8mb4_unicode_ci;
USE lims_db;
```
```

3.1 用户与权限管理 (users, roles, permissions)

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```
```sql
CREATE TABLE users (
 id INT PRIMARY KEY AUTO_INCREMENT,
 username VARCHAR(50) NOT NULL UNIQUE,
 password VARCHAR(255) NOT NULL,
 email VARCHAR(100) NOT NULL UNIQUE,
 full_name VARCHAR(100),
 role_id INT NOT NULL,
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE
CURRENT_TIMESTAMP,
 FOREIGN KEY (role_id) REFERENCES roles(id) ON DELETE CASCADE
);
```
```

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```
```sql
CREATE TABLE roles (
 id INT PRIMARY KEY AUTO_INCREMENT,
```

```
role_name VARCHAR(50) NOT NULL UNIQUE,
description TEXT
);
...
```

权限表 (permissions)

```
```sql  
CREATE TABLE permissions (  
    id INT PRIMARY KEY AUTO_INCREMENT,  
    permission_name VARCHAR(50) NOT NULL UNIQUE,  
    description TEXT  
);  
...
```

角色权限关联表 (role_permissions)

```
```sql  
CREATE TABLE role_permissions (
 role_id INT,
 permission_id INT,
 PRIMARY KEY (role_id, permission_id),
 FOREIGN KEY (role_id) REFERENCES roles(id) ON DELETE CASCADE,
 FOREIGN KEY (permission_id) REFERENCES permissions(id) ON DELETE CASCADE
);
...
```

## 3.2 样本管理 (samples)

样本表 (samples)

```
```sql
```

```
CREATE TABLE samples (  
    id INT PRIMARY KEY AUTO_INCREMENT,  
    sample_code VARCHAR(50) NOT NULL UNIQUE,  
    sample_name VARCHAR(100) NOT NULL,  
    sample_type VARCHAR(50) NOT NULL,  
    collected_by INT NOT NULL,  
    received_by INT NOT NULL,  
    received_date DATETIME NOT NULL,  
    storage_location VARCHAR(100),  
    status ENUM('待检测', '检测中', '已完成', '已销毁') DEFAULT '待检测',  
    created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,  
    updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE  
CURRENT_TIMESTAMP,  
    FOREIGN KEY (collected_by) REFERENCES users(id) ON DELETE SET NULL,  
    FOREIGN KEY (received_by) REFERENCES users(id) ON DELETE SET NULL  
);  
```
```

### 3.3 检测与实验管理 (experiments, experiment\_results)

实验表 (experiments)

```
```sql
```

```
CREATE TABLE experiments (  

```

```

id INT PRIMARY KEY AUTO_INCREMENT,
experiment_code VARCHAR(50) NOT NULL UNIQUE,
sample_id INT NOT NULL,
assigned_to INT NOT NULL,
experiment_type VARCHAR(100),
start_date DATETIME,
end_date DATETIME,
status ENUM('待处理', '进行中', '已完成') DEFAULT '待处理',
remarks TEXT,
created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
updated_at  TIMESTAMP  DEFAULT  CURRENT_TIMESTAMP  ON  UPDATE
CURRENT_TIMESTAMP,
FOREIGN KEY (sample_id) REFERENCES samples(id) ON DELETE CASCADE,
FOREIGN KEY (assigned_to) REFERENCES users(id) ON DELETE SET NULL
);
...

```

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```

```sql
CREATE TABLE experiment_results (
 id INT PRIMARY KEY AUTO_INCREMENT,
 experiment_id INT NOT NULL,
 result_data TEXT NOT NULL,
 result_status ENUM('正常', '异常') DEFAULT '正常',
 recorded_by INT NOT NULL,
 recorded_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 FOREIGN KEY (experiment_id) REFERENCES experiments(id) ON DELETE
CASCADE,
 FOREIGN KEY (recorded_by) REFERENCES users(id) ON DELETE SET NULL

```

---

);

...

### 3.4 设备与仪器管理 (equipment)

设备表 (equipment)

```sql

```
CREATE TABLE equipment (  
    id INT PRIMARY KEY AUTO_INCREMENT,  
    equipment_code VARCHAR(50) NOT NULL UNIQUE,  
    equipment_name VARCHAR(100) NOT NULL,  
    model VARCHAR(50),  
    manufacturer VARCHAR(100),  
    purchase_date DATE,  
    last_calibration DATE,  
    next_calibration DATE,  
    status ENUM('正常', '维修', '报废') DEFAULT '正常',  
    location VARCHAR(100),  
    assigned_to INT,  
    created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,  
    updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE  
CURRENT_TIMESTAMP,  
    FOREIGN KEY (assigned_to) REFERENCES users(id) ON DELETE SET NULL  
);  
...
```

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3.6 质量控制 (quality_control)

质量控制记录表 (quality_control)

```
```sql
```

```
CREATE TABLE quality_control (
 id INT PRIMARY KEY AUTO_INCREMENT,
 sample_id INT NOT NULL,
 experiment_id INT NOT NULL,
 qc_type VARCHAR(50) NOT NULL,
 qc_status ENUM('合格', '不合格') DEFAULT '合格',
 reviewed_by INT NOT NULL,
 reviewed_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 FOREIGN KEY (sample_id) REFERENCES samples(id) ON DELETE CASCADE,
 FOREIGN KEY (experiment_id) REFERENCES experiments(id) ON DELETE
CASCADE,
 FOREIGN KEY (reviewed_by) REFERENCES users(id) ON DELETE SET NULL
);
...
```

## 3.7 数据统计与报表 (reports)

报告表 (reports)

```
```sql
```

```
CREATE TABLE reports (  
    id INT PRIMARY KEY AUTO_INCREMENT,  
    report_code VARCHAR(50) NOT NULL UNIQUE,  
    experiment_id INT NOT NULL,  
    generated_by INT NOT NULL,  
    report_content TEXT NOT NULL,  
    created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,  
    FOREIGN KEY (experiment_id) REFERENCES experiments(id) ON DELETE
```

```
CASCADE,  
    FOREIGN KEY (generated_by) REFERENCES users(id) ON DELETE SET NULL  
);  
...
```

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